

Fluorescence Correlation Spectroscopy (FCS)

Homework TD


A- Scientific paper Review.


- 1- Read Material & Method of Guiziou S. et al. **Nucl. Acids Res.**(2016) doi: 10.1093/nar/gkw624.
- 2- Read the entire paper Zuniga A. et al. **ACS Synth. Biol.** 2021, 10, 3527–3536.
to understand the goal of N&B experiments on ALPAGA living cell in different growth conditions (+/- Glucose and +/- Lactate.)

- 1- Write precise protocol for sample preparation (growth media, OD, Immobilization,...)
- 2- Write data acquisition protocol (laser wavelength, power, size image, dwell time, # image,...)


A- Data analysis simulation with PaTrack software:

- Open PaTrack-better
- Load Test tif files (m_04-920 nm_6mW_40 usec_50 frames_Ch2
- Follow the instructions.






- Autodetect then delete inconsistent.
- Play with the threshold (cursor above Autodetect) to detect most of the bacteria
- Right click to delete, separate, ect cells....
- To add cell select the red bacteria on top then trace a line all along the bacteria that you want add using the mouse.



- Detect all cell to get as follow



Once you're done with cell selection, save your data.

A .spe file will be saved automatically in the same directory.

Open .spe file in bloc-notes

02-920 nm_6mW_40 usec_50 frames_Ch2.spe - Bloc-notes

Fichier Edition Format Affichage Aide

C:\Users\CClerte\Desktop\Rho Polymerase project_Manu_bossi\2015\04142015\strain 191 Turquoise Rho\02-920 nm_6mW_40 usec_50 fr

Channel	ID	FlagROI	X(px)	Y(px)	Ang.(°)	a(um)	b(um)	Area.Elli.(um2)	Area.Bact.(um2)	Vol.Elli.(um3)	Vol.Bact.(um3)
1	1	1	185.78	18.41	-98.46	15.60	4.60	225.46	268.91	1382.58	1870.12
1	2	1	175.33	18.09	-99.91	12.59	3.36	133.02	159.66	596.51	815.09
1	3	1	50.67	20.96	-158.39	21.54	3.84	259.72	318.05	1328.81	1874.88
1	4	1	204.57	23.86	-136.87	16.67	4.11	215.10	259.39	1178.25	1622.16
1	5	1	168.08	29.96	-99.76	17.27	3.57	193.70	235.67	922.24	1287.99
1	6	1	116.71	32.76	-93.06	17.57	3.77	208.26	252.95	1047.51	1458.83
1	7	1	86.68	28.01	-20.49	15.78	4.16	206.25	247.74	1144.31	1565.57
1	8	1	149.76	32.83	-39.74	15.06	3.77	178.24	214.75	895.45	1231.13
1	9	1	234.60	36.26	-134.60	14.81	3.89	181.04	217.52	939.00	1285.22
1	10	1	208.33	40.93	-135.45	20.26	3.71	236.32	289.06	1169.69	1647.40
1	11	1	90.19	42.35	-43.08	17.70	4.14	230.32	278.53	1271.71	1758.84
1	12	1	45.28	39.06	-160.67	19.36	4.09	248.94	302.59	1358.57	1894.24
1	13	1	191.10	41.85	-154.13	17.73	3.90	217.40	263.73	1131.42	1572.58
1	14	1	125.60	49.00	-52.52	13.56	3.86	164.18	196.29	843.92	1145.89
1	15	1	231.75	55.20	-119.52	20.14	3.69	233.61	285.74	1149.98	1619.57
1	16	1	96.75	55.28	-74.68	13.45	3.69	156.11	187.05	768.70	1047.55
1	17	1	184.35	53.46	-123.65	12.99	3.76	153.51	183.31	770.17	1043.68
1	18	1	149.40	56.05	-36.85	21.06	3.35	221.77	272.72	991.19	1407.90
1	19	1	52.91	53.07	-5.97	20.86	4.03	263.90	322.09	1417.23	1988.99
1	20	1	156.08	62.01	-35.87	23.93	3.33	250.36	309.25	1111.90	1590.45
1	21	1	211.10	58.81	-163.82	21.38	3.97	266.37	325.65	1408.78	1982.46
1	22	1	163.10	67.96	-35.64	20.07	3.57	225.28	275.88	1072.97	1513.98
1	23	1	62.04	70.49	-9.91	12.57	3.71	146.43	174.63	724.13	979.33
1	24	1	214.55	78.24	-68.46	12.57	3.48	137.49	164.65	638.28	869.02
1	25	1	136.01	82.42	-123.00	17.15	3.80	204.79	248.35	1038.03	1441.99
1	26	1	22.61	83.31	-100.49	14.72	3.66	169.27	204.02	826.27	1136.64
1	27	1	37.52	72.80	-4.04	11.61	3.73	135.85	161.06	674.81	903.93
1	28	1	83.85	80.98	-55.43	13.85	3.46	150.58	181.45	694.66	955.24
1	29	1	223.55	81.91	-72.79	11.76	3.58	132.30	157.44	631.69	851.35
1	30	1	127.70	88.83	-117.96	18.53	3.70	215.13	262.18	1060.25	1484.60
1	31	1	201.58	88.25	-78.26	15.33	3.64	175.40	211.94	851.79	1176.50
1	32	1	110.71	91.43	-105.86	17.07	3.58	191.96	233.42	915.91	1277.88

Ln 1, Col 1 100% Windows (CRLF) ANSI

Copy and paste n, e, B and <F> data in excel

- Then proceed with file Rho_02-920 nm_6mW_40 usec_50 frames_Ch2) in PaTrack.
- Calculate the average true brightness e (or ε) for both set of data.
- Discuss and Comment the results obtained.
- If Fluorescence intensity background is = 1; Calculate the corrected brightness and the corrected Number of molecule/cell.
- Think about and Write a script to calculate the absolute Number of Molecule/cell from N&B data that you'll collect during the practical.

B- Diffusion Coefficient Calculus Exercises.

- Calculate theoretical Diffusion Coefficient for Alexa488 and GFP in Water.
- Look for in biblio a experimental measured Diffusion Coefficient of EGF in Cytoplasm (either Eukaryote or Prokaryote). Calculate the apparent crowdedness of the cytoplasm.
- Write autocorrelation function as function of W_0 , Z_0 , τ and τ_D . What happens when $\tau = \tau_D$
- Calculate the diffusion Coefficient for a monomeric and a dimeric form of a protein (MW= 100KDa). If you have a perfect mixture of monomer and dimer in your tube, can you observe and quantify each species in a FCS autocorrelation experiment?

C- Simulations Point FCS data using SimFCS software.

Exercise 1: Simulation of single species FCS.

1-1: Setting up the problem

1_Open SimFCS (Fig. 1) (As administrator)




Fig 1

2_Select “Simulate FCS data” in the action panel at the left-hand side of the screen (Fig 2).




Fig 2

A “simform” panel opens up (Fig. 3).




Fig. 3

3_Select “Parameters” from the menu at the top of the screen (Fig. 3).

PSF: This tab has the options for the shape of the point spread functions (Fig. 5)




Fig 5

4_ Select “3-D Gaussian” from the pull-down menu.

5_ Set the **Radial Waist** (r or y) to **0.3** (the units are microns) and the **Axial Waist** to **1.5** μm . Ignore the value of **Radial Waist (x)**—this refers only to point spread functions that are not cylindrically symmetric should be set to zero.

The **plotting threshold** is simply for visualization and will not affect the outcome of the simulation.

6_ Set the **box size** to **64**.

This is actually the half size of the box. Since the simulation moves the particles in a random walk along a grid of points spaced $0.05 \mu\text{m}$ apart, this will give a box with $6.4 \mu\text{m}$ along each dimension.

7_ Set the **multiple points** field to **1**.

Ignore the **slit separation** as this is for more complex simulations, as well as **x,y,z for nonscanning**

and leave **Perfect 0_9N.psf**

Note that if you click the last line of this tab gives you the volume of the box and the size of the PSF in μm^3

Particles: This tab allows you to choose the numbers and types of particles for the simulation (Fig 6).

We will be doing a single species simulation, so only the first row will have particles.




Fig 6

8_ In column “# mol” column choose **100** molecules.

Note: The maximal number of particles is 4000.

9_ In column “Br1 (cpsm)” choose a brightness of **100,000** counts per second per molecule (cpsm).

This is the average intensity at the peak of the point spread function.

10_ In column “D ($\mu\text{m}^2/\text{s}$)” set the diffusion coefficient to **40** $\mu\text{m}^2/\text{s}$.

Note that at the bottom left hand corner of the tab, there is a Max D value in $\mu\text{m}^2/\text{s}$. This maximum value is set by the simulation based on the time per channel (clock frequency) which we will set later. This D value is the value at which the probability of making two grid steps in one-time step becomes significant. It is possible to set the diffusion coefficient to a value higher than this, but the simulation will be physically unrealistic, because grid steps greater than one are not allowed.

Ignore the **rotation value** for now.

11_ Set “Br2” to zero.

This is the brightness of the particle in the second channel. It is set to zero because we will only be simulating one detection channel.

The “**Flow**” column introduces a directional movement and will not be used in this exercise.

12_ Select type “1” for normal particles in the seventh column “1..6”.

At the bottom right of the table is a legend showing the type of particle that this last column can designate.

A “2” in this column indicates that the particle will rotate with a time constant equal to the value in the rotation column, which is in units of time channels. A “3” and “4” in this column refer to different shapes of macromolecular assemblies, like molecules on a stick (DNA) and molecules on a disk (rafts).

Misc: This tab contains options for rotation, flow, background, and particle confinement (Fig. 7).




Fig 7

13_ Select “randomize data” and

14_ “Use Poisson dist for counts.”

These settings are appropriate for simulating photon counting, not analog detection. For analog detection, on top of a Poisson distribution there is an additional exponential distribution due to the process of electron multiplication in the photomultiplier.

Time/scan: This tab allows for the simulation of circular scanning, line scanning, camera scanning or raster scanning FCS (Fig. 8).




Fig. 8

15_ Make sure the “**Scanning**” check box is **unchecked**.

16_ Set the clock **frequency** to **100,000 counts/sec/particle**. Note that in this simulation we are using 10 μ s dwell time.

The frequency is the inverse of the dwell time in seconds

17_ Set the “**Cycles**” field to **50**.

One cycle is the amount of time required to collect 2¹⁵ data points.

*At the right-hand side of the tab, the actual simulation time (**time(s)=**) is shown in seconds.*

The “**Binning Mode**,” “**Use sections**,” and “**Number of bins for linear scale**” fields are for display during the simulation and will not affect the outcome.

The following tabs will be ignored for now:

Stick/Disk: This tab allows the simulator to change the shape and properties of the stick or disk particles selected under the Particles tab.

Reaction: This tab allows for reactions (exchange of properties) between particles (rows) defined in the

Particles section. This allows to simulate bleaching, FRET or other chemical reactions.

Display: The tab allows for the change of the display during the simulation run. It does not affect the simulation.

TZone 2D: The tab allows for parameters related to transient confinement of particles.

Binding 2D: Used to introduce sites of binding in the simulation

Channel 2D: Used to simulate diffusion/transport across a channel

3D object: Used for the simulation of 3D shapes.

1-2: Running the Simulation:

18_ Click on “**Parameter**” (Fig. 9)




Fig 9

19_ Click on “**Start_simulation**” at the top of the simform window.

20_ Select “**Isotropic diffusion**” from the pull-down menu.

After a moment, the simulation will start.

If you like, select “Show PSF” (this will show a 3D threshold plot of the PSF). This can also be done before the simulation to demonstrate the different PSFs. In addition, the particles can be plotted in real time (“Plot particles” checkbox). All of these display options slow the simulation considerably, so it is advisable to turn them off for the majority of the simulation.

Most plots in SimFCS can be modified by double clicking on the graph. The plotting routine is TeeChart Pro and a free editor for Tee Chart Pro, Tee Chart Office, can be freely downloaded at:

http://www.steema.com/download/other_projects.

Also, during the simulation, the data will register on the main SimFCS screen (Fig. 9) similar to the way real data acquisition would register. The top plot, an accordion plot, displays the intensity per second in red, and the $G(0)$ value of channel 1, estimated by the variance divided by the intensity squared, in green. On the bottom plot, the logarithmically binned autocorrelation function is plotted. This function is calculated as the average of each cycle (215 points) and is calculated using the multiple tau approach. Note that this should only be used for visual purposes during the simulation. After the simulation is complete, the “Recalculate” button (on the lower right of the screen) can be pushed to get a more robust estimation of the autocorrelation function.

1-3: Analyzing the Data:

21_ Click “**Analysis**” under the main SimFCS screen, on the top menu (Fig. 10).




Fig 10

22_ Select “**Large vector correlation**” from the pull-down menu.
A new window (**corrform**) opens (Fig. 11).




Fig 11

This section of the program calculates the autocorrelation function based on the **fast Fourier transform (FFT)**. Since the FFT requires data to have a size that is a power of two, the autocorrelation function is calculated in appropriately sized sections. In this example:

23_ Select the entry **256K** (not 256!).

24_ Select “**Plot each segment**” (screen bottom left), before starting the calculation to view these segments as they are calculated.

25_ Press the “**calculate**” button to begin the calculations.

The individual autocorrelation segments can be deleted. This capability is extremely useful in deleting, for example, a section of the data that contains a spike, often due to a diffusing aggregate.

To delete an anomalous autocorrelation segment:

-Click on it when it will appear red.

-Press “Delete segment”, then

-Press “Recalculate average” to recalculate the autocorrelation average without the (offending) segment

Note that the errors are also calculated and can be visualized by selecting “Plot errors.” When the data are saved (“Save correlation files”), the errors are saved as well. They can then be used as weighting factors for fitting.

At the end of the calculation, the average is shown on the top graph (may be difficult to see if you are plotting each segment). The bottom graph shows the logarithmically binned data. This format is generally

accepted for FCS data. On the right-hand side of the page is the photon count histogram (PCH) and its average.

26_ Select the “fit” tab on top of the corrform window (Fig.12).




Fig 12

The “corrfit” window will pop up (Fig. 13).




Fig 13

27_ Input the values for the **waist=0.3** (in micron)

28_ Deselect 2-photon formulas.

29_ Initialize the fit by extrapolating the value of **G(0)** from the graph.

This step is not necessary, but it is important that the initial value of $G(0)$ is different from zero.

30_ Select “**use weights**”

31_ Set “**Fix Bkgd**” **unchecked**.

This will allow the background to vary.

32_ Click on “**perform fit**”.

QUESTION form exercise 1:

- Calculate the concentration of particle in the simulation box ($1\mu\text{m}^3 = 1\text{fL} = 10^{-15}\text{L}$)
- If you change beam waist or any other simulation parameter (diffusion, concentration, brightness, sample frequency). Comment your observations.
- Fit autocorrelation function obtained in the simulation to different radial waist. Comment results obtained.
- Fit autocorrelation function to different waist ratios (w_0/w_z). Comment on the sensitivity of FCS to the z dimension.

Exercise 2: Simulation of 2 species.

2-1: Setting up the problem:

For the formation of a dimer, the brightness and volume might both be expected to double. The diffusion coefficient is inversely proportional to the cube root of the particle volume (the particle radius), so it will be multiplied by a factor of ($0.5^{1/3} = 0.794$), or approximately $32\mu\text{m}^2/\text{s}$, if $40\mu\text{m}^2/\text{s}$ is

the diffusion coefficient of the monomer.

Set up the simulation as we did in exercise 1 (steps 1 to 7), except for the particles section.

1_ Open **SimFCS** (Fig. 1)

2_ Select “**Simulate FCS data**” in the action panel at the left hand side of the screen (Fig 2).

3_ Select “**Parameters**” from the menu at the top of the screen (Fig. 3).

4_ Select “**3-D Gaussian**” from the pull down menu (Fig. 5).

5_ Set the **Radial Waist** (r or y) to **0.3** (the units are microns) and the **Axial Waist** to **1.5 μm** (Fig. 5) .

6_ Set the **box size** to **64** (Fig. 5) .

7_ Set the **multiple points** field to **1** (Fig. 5) .

In the particles section (Fig. 15):

35_ Set the **number of particles** at **75** (mol#) with a **brightness** of **100,000** cpsm (Br1 cpsm)

36_ Set **25** particles at **200,000** cpsm.

37_ Set the **diffusion** at **20** and **100** (D ($\mu\text{m}^2/\text{s}$))




Fig. 15

2-2: Running the Simulation:

Run the simulation following Steps 18 to 20 of exercise 1 (Fig. 9).

18_ Click on **"Parameter"**.

19_ Click on **"Start_simulation"** at the top of the simform window.

20_ Select **"Isotropic diffusion"** from the pull-down menu

2-3: Analyzing the Data (SimFCS):

Repeat steps 21-32 of exercise 1 (Figs. 10 to 13)

21_ Click **"Analysis"** under the main SimFCS screen, on the top menu (Fig. 10).

22_ Select **"Large vector correlation"** from the pull-down menu (Fig. 10).

23_ Select **256K** (not 256!) (Fig. 11).

24_ Select **"Plot each segment"** (screen bottom left), before starting the calculation to view these segments as they are calculated (Fig. 11).

25_ Press the **"calculate"** button to begin the calculations (Fig. 11).

26_ Select the **"fit"** tab on top of the corrfom window (Fig.12).

27_ Input the values for the **waist=0.3** (Fig. 13)

28_ **Deselect 2-photon formulas.** (Fig.13)

29_ **Initialize** the fit by extrapolating the value of **G(0) from the graph** (Fig. 13).

30_ Select **"use weights"** (Fig.13).

31_ Set **"Fix Bkgd"** unchecked (Fig.13).

32_ Click on **"perform fit"** (Fig.13).

To fit the PCH data:

38_ Select on the upper left corner on the down arrow menu and switch from **"autocorrelation 1"** to **"PCH 1"**.

39_ Initialize the values for N and epsilon and leave N1, epsilon 1 and N2, epsilon 2 unchecked to allow the two species to vary.




Fig. 16

QUESTIONS FOR EXERCICES 2:

- Calculate the concentration and diffusion of each species.
- If you change beam waist and or any others parameters (diffusion Coefficients-concentration-Brightness-Sample frequency,...) Comment your observations.

Exercise 3. Simulation of two-color experiments: cross-correlation analysis of heterodimer formation

3-1: Setting up the problem:

40_Set up the simulation with three types of particles as seen in Fig 17 (arrows).




Fig. 18

25_ Press the “**calculate**” button to begin the calculations (Fig. 11).

26_ Select the “**fit**” tab on top of the corrform window (Fig.12).

27_ Input the values for the **waist=0.3** (Fig. 13)

28_ **Deselect 2-photon formulas.** (Fig.13)

29_ **Initialize** the fit by extrapolating the value of **G(0) from the graph** (Fig. 13).

30_ Select “**use weights**” (Fig.13).

31_ Set “**Fix Bkgd**” **unchecked** (Fig.13).

14th LFD Workshop in Advanced Fluorescence Imaging and Dynamics 24

32_ Click on “**perform fit**” (Fig.13).

To fit the PCH data:

33_ Click on the upper left corner on the down arrow menu and switch from “autocorrelation 1” to “PCH 1” (Fig.14).

34_ Click on “**perform fit**” (Fig.14).

Results are shown in Fig. 19.




Fig. 19

Autocorrelation in Channel 2.

Select "Ch2 autocorrelation" and do that analysis the same way as Ch1.

Note: Do not forget to properly **initialize** the fit by extrapolating the value of $G(0)$ from the graph.




Fig. 20

The autocorrelation for channel 2 should appear in blue in the bottom graph (Fig. 20). The information in the upper graph is overwritten when a new analysis is done, so make sure that you delete any outliers and recalculate the average before analyzing a new channel.

Cross-correlation between channel 1 and 2

Now select "Ch1-Ch2 cross-corr" and perform the analysis.

Note: Do not forget to properly **initialize** the fit by extrapolating the value of $G(0)$ from the graph.




Fig. 21

The cross-correlation will appear in green along with the two autocorrelations (Fig. 21). You can also do the "Ch2-Ch1 cross-corr" if you like, but the result will be the same and the Ch1 Ch2 data will be overwritten. Also, only one cross correlation can be saved at a time. The cross-correlation is only sensitive to species that diffuse together in both channels, so you can fit the data to a Gaussian single species model as in Exercise 1.

QUESTIONS FOR EXERCICES 3:

- If you change beam waist and or any others parameters (diffusion Coefficients-concentration-Brightness-Sample frequency,...) Comment your observations.

Equations:

$$V_{\text{PSF}} = \int \text{PSF} d\mathbf{r}^3$$

$$\gamma = \frac{\int \text{PSF}^2(\mathbf{r}) d\mathbf{r}^3}{\int \text{PSF}(\mathbf{r}) d\mathbf{r}^3}$$

3D Gaussian:

$$\text{PSF}(r, z) = \exp\left(\frac{-2r^2}{\omega_x^2} + \frac{-2z^2}{\omega_z^2}\right)$$

$$V_{\text{3DG}} = \omega_x^2 \cdot \omega_z \cdot \left(\frac{\pi}{2}\right)^{3/2} = 1.9687 \cdot \omega_x^2 \cdot \omega_z$$

$$\gamma_{\text{3DG}} = 2^{3/2}/8 = 0.3536$$

$$G(0) = \sum_m \left(\frac{I_m}{\langle I \rangle} \right)^2 G(0)_m, \text{ where } I_m \text{ is the intensity of species } m.$$

$$\langle I_m \rangle = \epsilon_m \cdot N_m, \text{ where } \epsilon_m \text{ and } N_m \text{ are the brightness and number of species, } m.$$

$$\langle I \rangle = \sum_m \epsilon_m \cdot N_m.$$

Equations:

$$G_m(\tau) = \frac{\langle I_{a1}(t) \cdot I_{a2}(t+\tau) \rangle}{\langle I_{a1}(t) \rangle \langle I_{a2}(t) \rangle} = \frac{\langle I_{a1}(t) \cdot I_{a2}(t+\tau) \rangle}{\langle I_{a1}(t) \rangle \langle I_{a2}(t) \rangle} - 1 \approx \frac{\langle I_{a2}(t) \cdot I_{a1}(t+\tau) \rangle}{\langle I_{a1}(t) \rangle \langle I_{a2}(t) \rangle} - 1.$$

where $ch1$ and $ch2$ refer to channel 1 and channel 2 respectively. $G_m(0) = \gamma \sum_{m=1}^M \frac{z_{m,ch1} \cdot z_{m,ch2} \cdot N_m}{\langle I_{ch1} \rangle \langle I_{ch2} \rangle}$,

where N_m is the number of type m particles in the focal volume and $z_{m,ch1}$ and $z_{m,ch2}$ are the brightnesses of species m in each channel.

$$\langle I_{ch1} \rangle = \sum_m z_{m,ch1} \cdot N_m.$$